

Performance of the OMEGA EP High-Energy Short-Pulse Laser System



B. E. Kruschwitz *et al.*
University of Rochester
Laboratory for Laser Energetics

**International Conference on
Ultrahigh Intensity Lasers
Watkins Glen, NY
26 September – 1 October 2010**

Summary

OMEGA EP is operational and advancing toward its rated performance level



- **OMEGA EP is a productive high-energy short-pulse laser enabling scientific advances**
- **LLE is working to extend the performance to full specifications**
 - **short-pulse optics performance is being improved—
knowledge of safe operating limits is being developed**
 - **apodization is evolving to optimize beam uniformity—
active beam correction is being implemented**
- **Advancements in short-pulse diagnostics have resulted in improved on-shot characterization capabilities and enhanced system performance**


Collaborators




**C. Dorrer, M. J. Guardalben, J. H. Kelly, M. D. Moore, J. Qiao, L. J. Waxer,
I. A. Begishev, J. Bromage, S.-W. Bahk, L. Folsbee, D. Irwin, S. D. Jacobs,
R. Jungquist, T. J. Kessler, R. W. Kidder, S. J. Loucks, J. R. Marciante,
R. L. McCrory, D. D. Meyerhofer, S. F. B. Morse, A. V. Okishev, J. B. Oliver,
G. Pien, J. Puth, A. L. Rigatti, T. C. Sangster, A. W. Schmid, M. J. Shoup III,
C. Stoeckl, K. A. Thorp, and J. D. Zuegel**

**University of Rochester
Laboratory for Laser Energetics**

The Extended-Performance (EP) addition to OMEGA is located adjacent to the OMEGA Laser System




OMEGA EP is an exceptionally flexible high-energy, high-power laser facility




- OMEGA EP delivers two separate kilojoule-level, picosecond-pulse beamlines to the OMEGA EP target chamber
- The two short-pulse beams will be capable of co-propagating to either the OMEGA or OMEGA EP target chamber
- OMEGA EP delivers nanosecond UV pulses in four beamlines to the OMEGA EP target chamber
- The kilojoule-level, nanosecond UV beams can be used together with the short pulse beams

OMEGA EP has performed >1000 target shots for a variety of experiments

- Target shots to date
 - 750 short pulse only (188 to Ω TC)
 - 200 long pulse only
 - 78 combined short and long pulse
- Scientific advances enabled
 - integrated fast ignition
 - positron production
 - **proton radiography**





Proton radiography measurements of target charging




P. M. Nilson *et al.*, "Scaling Hot-Electron Generation to High-Power, Kilojoule-Class Laser-Solid Interactions," submitted to Phys. Rev. Lett.

OMEGA EP IR beamlines use a folded architecture based on the NIF




Each compressor is comprised of four tiled-grating assemblies (TGA's)

Tiled-grating assembly




Short-pulse diagnostic path



- The TGA's are interferometrically tiled inside the grating compressor chamber at vacuum
- Fourier fringe analysis is used to tile the TGA's and to retrieve TGA wavefront
- Differential piston, tip, and tilt are automatically calculated and removed for initial tiling
- Tiled positions are maintained by closing the actuator control-loop with feedback of position displacement sensors
- Submicroradian angular stability is achieved


An energy ramp to >2.1 kJ at ~12 ps was performed in September 2009




Damage was observed on the final grating and downstream optics. Energy is currently limited to 1 kJ at ≥ 10 ps as a result.

Typical performance of a BL2 1-kJ shot (need to come up with a sentence title)

Near field




Spectrum




Pointing $\sim 15\text{-}\mu\text{m}$ rms
Energy contrast $\sim 10^4$
Power contrast $\sim 10^6$
for 10-ps pulses

Wavefront



Target-plane focal spot



OMEGA EP is aggressively working to meet its design specifications



Design Energy Performance	Current Energy Performance
1 PW (600 J at 0.6 ps)	0.5 PW (300 J at 0.6 ps)
2.6 kJ for $\tau \geq 10$ ps	1.0 kJ for $\tau \geq 10$ ps

Challenges

- Increase short-pulse damage threshold of gratings and transport optics
- Optimize beam fill factor and minimize beam modulation
- Minimize on-target focus-spot size
- Minimize prepulse


LLE has a comprehensive program to improve the performance of short-pulse optics



- Dedicated short-pulse vacuum damage-testing facility
- *In-situ* grating damage detection
- MLD coating facility with electron-beam **or** plasma-assist deposition on large OMEGA EP size optics (~0.9-m diagonal), which provides options for environmental stability and/or film stress control
- Metrology
 - AFM and Nomarski microscopes that allow evaluation of full-scale OMEGA EP optics
 - high-resolution SEM for witness-sample assessment
- Theoretical modeling of electric-field enhancement to guide grating design
- Analytical support for cleaning chemistry studies




LLE built a dedicated vacuum, short-pulse-laser damage-test facility



- 1.053- μm ultrafast laser system provides tens of millijoules at a 0.1-Hz repetition rate with a 0.6- to 100-ps pulse-width range
- Vacuum vessel provides damage testing at 10^{-6} Torr

OMEGA EP uses a grating inspection system (GIS) to detect damage on the final grating assembly

Grating inspection system



Scanned image of TGA4





Contrast enhanced to emphasize damage features

- A line-shape illumination generator projects a focused line onto the grating surface
- A 10-bit linear charge-coupled device (CCD) measures the scattered light from any features on the grating surface
- The line-shape generator and the CCD are scanned across the grating to generate a composite two-dimensional image of the entire grating surface
- Features $\geq 250 \mu\text{m}$ in size are detected

OMEGA EP is evolving its capability to compensate gain nonuniformity and optimize beam quality


Current approach: static beam shaping with an apodizer




Future enhancement: a programmable spatial-light modulator will refine the statically shaped beam






Hamamatsu LCOS-SLM
Area: $12 \times 16 \text{ mm}^2$
Resolution: 600×792 ($20\text{-}\mu\text{m}$ pixel)
Dynamic range: 2 waves




The short-pulse diagnostics package (SPDP) is used for pre-shot and on-shot pulse characterization




- Temporal diagnostics
- Spectrometer
- Near field
- Far field
- Wavefront
- Energy
- On-target focal spot
- Contrast



The focal-spot diagnostic (FSD) uses remote, on-shot near-field characterization to determine the on-target energy density



A two-step strategy using phase retrieval from far-field image data is used to improve measurement accuracy



Step 1:

Produce accurate prediction of the focal spot measured by the nearby far-field CCD camera


1. Phase retrieval of static lens errors
2. Phase retrieval of average relative phase among beam segments
3. Fourier analysis to estimate effects of polychromatism

Step 2:


Obtain accurate focal-spot measurements in the distant target chambers

1. Phase retrieval of static calibration errors


A modal phase-retrieval algorithm is used to retrieve wavefront information from focal-spot images



FSD accuracy was confirmed with direct target-plane microscope images in low-energy qualification testing




Statistics for all low-energy FSM shots



OMEGA EP currently uses two diagnostics to measure temporal contrast


- Two operational diagnostics
 - scanning cross-correlator for 5-Hz OPCPA front-end characterization (Sequoia, Amplitude Technologies)
 - fast photodiodes and oscilloscope for on-shot nanosecond contrast characterization of BL1 and BL2



- One diagnostic under development: on-shot cross-correlator for contrast measurement in the 500-ps window before the main pulse

The contrast diagnostics have contributed to improved pulse contrast on both beamlines

Contrast of OPCPA at best compression



OMEGA EP is operational and advancing toward its rated performance level



- **OMEGA EP is a productive high-energy short-pulse laser enabling scientific advances**
- **LLE is working to extend the performance to full specifications**
 - **short-pulse optics performance is being improved—
knowledge of safe operating limits is being developed**
 - **apodization is evolving to optimize beam uniformity—
active beam correction is being implemented**
- **Advancements in short-pulse diagnostics have resulted in improved on-shot characterization capabilities and enhanced system performance**