

Demonstration of High-Contrast, High-Intensity Laser Pulses using an OPCPA Preamplifier in a Double CPA Ti:sapphire Laser System



Hiromitsu Kiriya

Laser Research and Development Group, Quantum Beam Science Directorate
Advanced Photon Research Center (APRC), Japan Atomic Energy Agency (JAEA)
8-1 Umamidai, Kizugawa, Kyoto 619-0215, Japan

International Conference on Ultrahigh Intensity Lasers (ICUIL)
Development, Science and Emerging Applications
October 28, 2008
Shanghai-Tongji, China

Contributors



Michiaki Mori, Yoshiki Nakai, Takuya Shimomura, Manabu Tanoue,
Hajime Okada, Shuji Kondo, Shuhei Kanazawa, Izuru Daito, Akito
Sagisaka, Jinglong Ma, Hiroyuki Daido, Sergei Bulanov, Toyoaki Kimura,
Toshiki Tajima

Kinkakuji-Temple
40 km



J-KAREN laser



Laboratory



Todaiji-Temple
2 km



Hawaii
7,000 km



Outline



- ✓ Objective
- ✓ “J-KAREN” high intensity laser at APRC, JAEA
- ✓ OPCPA preamplifier
- ✓ Cryogenic-cooled Ti:sapphire power amplifier
- ✓ Results on amplification, compression, and contrast performance
- ✓ Summary


Problem remains to be solved for the application of ultra-high intensity lasers in high-field physics

✓ Principal issue with multi-terawatt laser experiments

Modern Ti:sapphire chirped-pulse amplification (CPA) lasers reach **intensities greater than 10^{20} W/cm²**. However, in the laser systems, a background of the amplified spontaneous emission (ASE) can generate unwanted plasmas before the main pulse arrives on the target.

✓ Objective


To develop **$>10^{10}$ temporal contrast** laser supporting multi-terawatt power level with an OPCPA preamplifier in Ti:sapphire laser system



“J-KAREN” laser schematic




J-KAREN; JAEA Kansai Advanced Relativistic ENgineering



H. Kiriya *et al.*, Opt. Lett. 33 (2008) 645.
 H. Kiriya *et al.*, submitted to Opt. Commun.

Physics of parametric amplification

- ✓ Uses a monochromatic laser to provide broadband amplification



A. Dubietis *et al.*, *Opt. Commun.*, 88 (1992) 437.
I. N. Ross *et al.*, *Opt. Commun.*, 144 (1997) 125.

Spatial overlap

+

Synchronized pulses


+

Phase matching




Energy transfer

Our front-end can produce cleaned high energy ($> \mu\text{J}$) seed pulse for injecting into subsequent amplifiers



OPCPA technique is used for a high-contrast and broad bandwidth preamplifier


✓ OPCPA preamplifier scheme




Our OPCPA system amplifies the seed pulse to ~ 10 mJ with an amplification gain of over 10^4 and is shown to be well suited for broad band amplification



✓ Amplification gain




✓ Amplified spectrum




Signal pulse from OPCPA system is amplified to 280 mJ in 4-pass Ti:sapphire preamplifier

✓ View of Ti:sapphire preamplifier



✓ Extraction energy




We cool the Ti:sapphire crystal in power amplifier down to below 100 K, in order to increase its thermal conductivity while reducing the dn/dt, for negligible thermal focusing



✓ View of Ti:sapphire power amplifier



✓ Temperature of the crystal



✓ Thermal focal length


$$f_{th} = \frac{\kappa}{QL} \left(\frac{1}{2} \frac{dn}{dt} \right)^{-1}$$

Thermal conductivity
40 W/mK @ 300K
1000 W/mK @ 100K

dn/dt
1.3×10⁻⁵ @ 300K
2.5×10⁻⁶ @ 100K

Thermal focal length is measured to be over **3.8 km** at the maximum pumping condition (6.5 J, 10 Hz)

✓ Experimental setup



✓ Thermal focal length


$$f_{th} = d - f \left(1 - \frac{f}{\Delta S} \right)$$

Measured thermal focal length; > 3.8 km !!


Using a cryogenic-cooled Ti:sapphire power amplifier, we amplify further to **3.2 J** and obtain good power stability of **$\pm 2\%$** (36,000 consecutive shots)



✓ Extraction energy




✓ Output power stability



Measured pulse duration is 30 fs with the energy of 1.7 J, corresponding to the peak power of 60 TW and temporal contrast is better than 10^{-10} at 10 Hz




✓ Pulse duration



The peak power is 60 TW (30 fs/1.7 J)

✓ Third order cross correlation




The actual contrast is found to be $\sim 10^{-11}$ level

Our J-KAREN laser system generates focused peak intensities in excess of 10^{20} W/cm² at 10 Hz




✓ Near-field spatial profile



Homogeneous and uniform spatial distribution

✓ Far-field spatial profile



8 μ m x 6 μ m (at the 1/e² points)

Peak intensity; > 10^{20} W/cm²

Conclusions



✓ We got over 10^{10} temporal contrast with over 10^{20} W/cm² intensity

J-KAREN laser of today;

10^{10} temporal contrast

10^{20} W/cm² intensity

10 Hz repetition rate

60 TW peak power

✓ We were encouraged by very successful upgrade

✓ We are currently supplying to laser-plasma accelerator experiments

HAMAMATSU



USHIO




TOSHIBA

PENTAX




SHIMADZU




by Prof. P. Bolton, Oct. 31

J-KAREN laser of tomorrow ??



✓ Outstanding issues for supplying PW level pulses -




150-J IR glass laser



60mm
Frequency doubler
80-J green laser



80mm
Ti:sap. crystal



Compressor chamber

Thank you very much for listening !



Acknowledgements for many fruitful discussions and stimulating inputs !

✓ Japan Atomic Energy Agency

Laser Research and Development Group, APRC

M. Tanaka, Y. Ochi, H. Sasao, R. Tateno, K. Kondo, A. Sugiyama

Photo-Medical Research Center (PMRC)

P. Bolton, M. Suzuki

✓ Institute of Laser Engineering, Osaka University

N. Miyanaga, K. Sueda

✓ Central Laser Facility, Rutherford Appleton Laboratory

J. Collier, C. Hernandez, C. Hooker, K. Ertel

✓ Laboratori Nazionali di Frascati dell'INFN

F. Terranova, F. Pegoraro